



MASTER YOUR PULSE

with

Accurate Mini
the World's Smallest
Multi-Vital Sign Monitor

Class II FDA APPROVED



Master Your Pulse

We all know people in our families or friends who are troubled, or even taken, by heart diseases. In fact, according to the World Health Organization (WHO), cardiovascular diseases are responsible for 17.9 million, near 1/3 of global deaths in 2019, and 85% of the deaths are caused by heart attack and stroke.

Therefore, stay aware of our health by monitoring blood pressure and other vital signs actively is the key to have a healthier life. Get an Accurate Mini, and master your pulse in style!



The Ultra Portable

Accurate Mini

Smallest Cuffless Calibration-free Multi-Vital Sign Monitor

Specification

Classification: Class II Medical Device
Dimension (Approx.): H37 x L66 x D18 mm
Weight (Approx.): 36g (without lanyard)
Display: 1-inch Full color TFT LCD
Battery: Enovix 3D Silicon™ Lithium-ion cell,
340 mAh, non-removal
Charging: USB Type-C
Performance: 500 measurements per charge
Memory: 500 set of measurements

Accurate Meditech Inc.

📍 8F-10, No. 12, Ln. 609, Sec. 5, Chongxin Rd.,
Sanhong Dist., New Taipei City, 241406,
Taiwan.
☎ +886-2-2999-1596
✉ info@accurate-meditech.com
🌐 <https://accurate-meditech.com>

Class II FDA APPROVED

2017 - 2023 © Accurate Meditech Ltd. All Rights Reserved.

