



# INSTANT Blood Pressure Reading Anytime Anywhere

510K Approval Cuffless Technology in your pocket.  
Your health can now be tracked more easily than ever.



Calibration  
Free



Instant



Silent



Accurate



microCare

	<b>Accurate Meditech</b>	<b>Apple Watch</b>	<b>Product Type</b>
<b>Product Type</b>	Clip	Watch	Wristband
<b>Application</b>	“Clip-up” to Smart Watch and Health Band	Wrist worn	Wrist worn without Smart Watch Function
<b>Output Type</b>	Blood Pressure Reading(mmHg)	Notification Only (Alerts)	Blood Pressure Reading(mmHg)
<b>Measuring Mode</b>	Instant Measurement	Trend tracking	Periodic Measurement
<b>Monitoring and Follow up</b>	Yes	No	Yes
<b>FDA Status</b>	Est. 2026 Q4 (OTC)	510(k) Clearance (alert noly)	510(k) Clearance / CE
<b>Positioning</b>	Medical-grade/self-monitoring blood pressure machine	High-End	High-End Relative to Features
<b>Calibration</b>	Calibration Free	N/A	Every 30 days
<b>Method</b>	Piezo / NIRS	PPG	PPG
<b>Price Level</b>	Affordable	High-End	High-End Relative to Features
<b>Customization</b>	Yes	N/A	N/A
<b>Connection</b>	Yes, api available	No	No